

Seiko's back surface, Free-form lenses offer vastly superior optical benefits with unequalled levels of optical performance, thinness, lightness and comfort. They also rate higher for patient satisfaction, compared to conventional progressive designs, particularly because free-form technology allows the production of completely customized lenses that specifically conform to every highly individual patient prescription.

With an exclusive license from free-form pioneer, Seiko Optical, Plastic Plus was the first Lab in North America to introduce these revolutionary Free-form PALs. And in only a few short years since, as the technology continues to evolve, and new generations of products are brought to market, Plastic Plus continues to lead in industry "firsts".



SEIKO Superior (*Best*)

- The most customized lens for all lifestyles
- 3 design options (Near, Far, Balanced)
- 11 corridor lengths (8mm to 18mm)
- 51 options for near zone inset
- 21 options for pantoscopic tilt
- Advanced Aspheric Compensation improves clarity and visual comfort
- Multi-polar Astigmatic correction improves panoramic vision and image stability
- Add powers from +0.50 to +4.00 (8 & 9mm corridors to +3.00)



SEIKO Surmount (*Better*)

- SEIKO Surmount - Soft design with the best vision in all zones
- SEIKO Surmount WS - Hard design for wider viewing zones
- Convex add curve on back surface for slimmer plus Rx's
- Ideal for difficult Rx's with very high sphere, cylinder, adds, and or prism
- Variable inset based on distance Rx, PD, and customized near reading distance
- As-Worn technology allows for fully compensated measured power throughout
- Multi-polar Astigmatic correction improves panoramic vision and image stability



SEIKO Supernal (*Good*)

- A good choice for all wearers, combines soft and hard design technology
- Developed using the new 3D Virtual Reality System for realistic simulation evaluation
- Non-linear corridor expands clear vision around the fitting point by 50%
- Automatic semi-variable inset based on total Rx
- Advanced aspheric compensation throughout improves clarity and visual comfort
- Multi-polar Astigmatic correction improves panoramic vision and image stability



SEIKO Supercede (*Standard*)

- Universal blended design provides wide distance and extra wide reading area
- Perfect spherical curve eliminates front surface distortion
- Placing optics on the back closer to the eye expands intermediate and near vision
- Advanced aspheric compensation in the intermediate and near zones

SEIKO

The First, Best, & Future of Free-Form

Branded PAL Design Availability

Available Lens Materials

Material	Clear HC	Polarized*		Transitions [®]		Transitions [®] XTRActive (Gray)	Transitions [®] Vantage
		Gray	Brown	Gray	Brown		
<input type="checkbox"/> 1.74	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/> 1.67	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
<input type="checkbox"/> 1.60	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>		
<input type="checkbox"/> 1.59	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> 1.53	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> 1.50	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>

Plastic Plus Premium Anti-Reflective Coating Options

Blue/Violet Filter

- The Blue Violet Filter Coating reduces reflection while blocking out potentially harmful HEV light caused by LED devices such as computers, laptops, tablets and phones

PP-S Super Oleophobic


- A true state-of-the-art coating featuring more anti-static, anti-reflective and durability layers.
- Super-Oleophobic Topcoat ensures visual impediments such as dust, water, oils, and fingerprints won't stick, and makes cleaning so easy

PP-O Oleophobic

- The same anti-static and anti-reflective properties as the premium PP-S AR with a slightly less slippery Oleophobic Topcoat.
- Provides outstanding scratch protection at a very competitive price point

PLASTIC PLUS

www.plastic-plus.com



FOR MORE INFORMATION AND PRODUCT AVAILABILITY CONTACT US TODAY!

TOLL FREE

Ontario 1.800.268.1617

Quebec 1.800.263.9417

Canada 1.800.387.8539


Fax: 1.800.530.9172

Tel: 416.789.4307

Fax: 416.789.7143

www.plasticplus.ca

Production Range



Canadian-owned and operated
SERVING CANADA SINCE 1974